

EPS HEP 2015, 22. - 29. 7. 2015, Vienna



Particle Identification Devices at the Belle II experiment

Rok Pestotnik

Jožef Stefan Institute, Ljubljana

on behalf of the Belle II collaboration

Outline

Introduction

Time of Propagation Counter

Aerogel RICH

Summary





Belle II at SuperKEKB B factory

Precision measurements of rare B, D and τ decays

- Other Belle II contributions at the conference:
 - J.Wiechczynski, The Belle-II Experiment at the Super-KEKB Collider, 25/7 10:20@
 Flavour Physics
 - Chunhua LI, The data acquisition and trigger system of the Belle II experiment, 25/7 10:00 @ Detector R&D
 - G.Casarosa, Inner tracking devices at the Belle II experiment, 24/7 15:00@ Detector R&D
 - T. Schlüter, Track Fitting in Belle II: the GENFIT Library and its Performance, 25/7 10:30, Detector R&D
 - M. Valentan, The Belle II Pixel Detector: How to deal with high occupancy, Detector R&D poster
- SuperKEKB will deliver 40 times higher event rates than KEKB.
- Belle II will collect 50ab⁻¹ by 2022 → around 10¹⁰ B mesons
 - Higher rates
 - Much higher backgrounds



Significant upgrade of the Belle detector needed.




Highly efficient π/K separation for momentum range up to 4GeV/c is needed to identify particles from various decays



Particle identification in Belle II








Barrel PID: Time of Propagation counter (TOP)

Cherenkov photons emitted in the quartz radiator \rightarrow **internally reflected** \rightarrow registered at the end of the bar by a fast position sensitive detector of single photons.





 $K/\pi \rightarrow Different \theta_c \rightarrow Different path length \rightarrow Different Time Of Propagation$



 θ_{c} reconstructed from:

hit position (x,y) and Time Of Propagation of photon.










TOP Design



Radiator: Quartz bar of high quality → Cherenkov photons reflect more than 100x


Photodetector → Hamamatsu SL10

- 4x4ch 10 µm Micro Channel Plate PMT
- \rightarrow works in B=1.5T, gain 2x10 6 , σ ~40ps, Q.E. min 24% @ 380nm, coll. eff. 50-55%
- \rightarrow 2x16 sensors /module \rightarrow 512 total

400 µm MCP

Front End Electronics:

High-speed 4GHz waveform sampling ASICs ("IRSX") 8 channels/chip σ< 50ps



27.5 mm

27.5 mm






TOP – Beam test @ SPring-8 LEPS 2013

Proof of the principle

Full scale prototype (with backup electronics)

2GeV/c positron beam


Excellent data and MC overlap







TOP Challenges




MCP PMT lifetime:


Atomic Layer Deposition to protect the photocathode from the positive ions

→ 284 conventional (will be replaced after degradation) and 284 ALD MCP-PMTs produced

High quality quartz bar:

- Flatness < 6.3 μm
- Roughness < 0.5 nm (RMS)
- Perpendicularity < 20 arcsec S1⊥S3,S4
- Parallelism < 4 arcsec S1|| S2







Front end electronics:

to read 8192 MCP PMT ch.

- 8 channels/ASIC chip,
- 4 ASICs / Carrier board,
- 4 Carrier boards / 8 PMT










TOP Production – 5 modules finished








Endcap PID: Proximity focusing Aerogel RICH

Goals and constraints:

- > 4 σ K/ π separation @ 1-3.5 GeV/c
- limited available space ~280 mm
- radiation tolerance (n,γ)

Selected type:

proximity focusing aerogel RICH





- <n> ~ 1.05
- $\theta_c(\pi) \approx 307 \text{ mrad } @ 3.5 \text{ GeV/c}$
- $\theta_{c}(\pi)$ $\theta_{c}(K)$ = 30 mrad @ 3.5 GeV/c
 - pion threshold 0.44 GeV/c,
 - kaon threshold 1.54 GeV/c
- neutron fluence: up to ~10¹² n/cm2
- radiation dose: up to ~1000 Gy

^{*} to increase the number of photons without degrading the resolution







ARICH - Aerogel radiator


Highly transparent (transmission length @400nm >30mm) hydrophobic silica aerogel

- Two 2cm thick layers $n_1 = 1.045$ $n_2 = 1.055$
- large tiles 18x18x2 cm³ to minimize photon losses at the edges



- 449 tiles already produced
- ready for cutting by water jets














ARICH - 144 ch. Hybrid Avalanche Photo-Detector











Last update: 2015/06/19

- HAPD developed with the Hamamatsu Photonics
- Works in a magnetic field of 1.5T
- Excellent separation of single photoelectrons







~ 80%(+10%) OK + 10% will be replaced



ARICH - Readout electronics




trigger

lto back-end

FPGA


limited space (~5cm) behind the HAPD

Board with 4 ASICs and Spartan6 FPGA







ARICH Beamtest @ T4-H6 CERN SPS 2011




- 6 HAPDs arranged as in the final design
- front-end electronics ~ close to final

> 5.5 σ K / π separation at 4 GeV/c












Momentum [GeV/c²]

Expected performance

Detailed simulation performed in the Belle II software framework – BASF2

results from release-00-05-00: ccbar events with nominal background



0.1

R.Pestotnik, PID at Belle II @ EPS-HEP2015





Summary

For efficient particle identification at Belle II two RICH detectors will be installed

- Time of propagation counter → Barrel PID
- Proximity focusing RICH with an aerogel as a radiator → Forward Endcap PID

TOP:

- 1st TOP module is under test in the cosmic ray setup, 4 have optic elements already glued, the rest will be produced by May 2016
- Combined TOP + CDC test is under preparation

ARICH:

- Highly transparent hydrophobic aerogel tiles (n=1.045,1.055) have been produced
- 353 out of 475 Hybrid avalanche photo detectors have been produced and tested
- Mass production of the ARICH readout boards is planned for the next month

All the key components of both detectors have been extensively tested on the bench and in the test beams.

The installation in the Belle II spectrometer → summer 2016

The beam test results and the detector simulations show that the designed detectors will have an excellent kaon identification efficiency at a rather low pion misidentification probability.










Belle II spectrometer




Belle II Detector







Schedule








Particle identification requirements

- improve π/K separation for few-body B decays
- good π/K separation for

$$b \rightarrow d \gamma, b \rightarrow s \gamma$$


- improve purity in fully reconstructed B decays
- low momentum (<1GeV/c) $e/\mu/\pi$ separation (e.g. for B \rightarrow KII)
- keep high the efficiency for tagging kaons





Design of Aerogel RICH








Performance in the magnetic field

- distortion of electric field lines at HAPD edge produces irregular shapes of areas covered by each channel
- in magnetic field photoelectrons circulate along the magnetic field lines and distortion disappears



- photoelectrons are backscattered from APD and produce long range cross-talk
- in magnetic field photoelectrons follow magnetic field lines and fall back on the same place

(radius of 10 keV electron in 1.5 T is~100µm)

scattered photons remain

